

System Description

- Flat substrate, xyz stage, inkjet deposition system
- Low cost, user- llable piezo-based inkjet print cartridges
- Built-in drop jetting observation system
- Fiducial camera for substrate alignment and measurement
- Variable jetting resolution and pattern creation PC-controlled with Graphical User Interface (GUI) application software
- Capable of jetting a wide range of uids
- Heated vacuum platen
- Cartridge cleaning station
- Includes software

Mechanical System

- Printable area
 - Substrate < 0.5 mm thickness: 210 mm x 315 mm (8.27 in x 12.4 in)
 - Substrate 0.5 25 mm thickness: 210 mm x 260 mm (8.27 in x 10.2 in)
- Repeatability: ± 25 μm (± 0.001 in)
- Substrate holder
 - Vacuum platen
 - Temperature adjustable; ambient to 60° C
- System Footprint: 673 mm x 584 mm x 419 mm (26 in x 23 in x 16 in)
- Weight approximately 50.7 kg (111.8 lbs)
- Power 100-120/200-240 VAC 50/60 Hz 375 W maximum
- Operating range 15-40° C at 5-80% RH non-condensing
- Altitude up to 2000 m
- Safety and EMC compliance
 - Safety: NRTL Certi ed to EN 61010-1, UL 61010-1, CSA 22.2 No. 61010-1
 - EMC: EN61326-1 Class A, FCC Part 15 Class A

Fiducial Camera

- Allows substrate alignment using reference marks
- Allows positioning a print origin or reference point to match substrate placement
- Provides measurement of features and locations
- Provides inspection and image capture of printed pattern or drops
- Provides cartridge alignment when using multiple cartridges
- Allows matching drop placement to previously patterned substrate

Cartridge

- Type: Piezo-driven jetting device with integrated reservoir and heater
- Usable Ink Capacity: Up to 1.5 ml (user- llable)
- Materials Compatibility: Many water-based, solvent, acidic or basic uids
- Not Available after August 30, 2019
 - Nozzles: 16 nozzles, single row, 100 dpi
 - Drop Volume: 1 (DMC-11601) and 10 (DMC-11610) picoliter nominal
- New Cartridge Products to be available in 2020. Contact FUJIFILM Dimatix for more information

Control PC and Application Software

- Pre-loaded patterned templates
- Pattern preview
- Editors: Pattern, piezo drive waveform, cleaning cycle, substrate setting
- Bitmap (1 bit) les accepted

Replaceable Items

- Print cartridge with one-time user- llable reservoir
- Cleaning station nozzle blotting pad
- Drop watcher uid absorbing pad

Dimatix Materials Printer DMP-2850





Drop Watcher View



Fiducial Camera View





www.dimatix.com

PDS00085 Rev.05 3/11/19



The Dimatix Materials Printer (DMP) is a benchtop materials deposition system designed for micro-precision jetting a variety of functional fiuids onto virtually any surface, including plastic, glass, ceramics, and silicon, as well as fiexible substrates from membranes, gels, and thin films to paper products. A complete turnkey system, the DMP facilitates developing and testing manufacturing processes and product prototypes. It also can be used for prototyping of products from fiexible circuits, RFID tags and displays to DNA arrays, and wearable electronics. By employing inexpensive exchangeable cartridges that researchers can fill with their own fiuid materials, the DMP system minimizes waste of expensive fiuid materials, thereby eliminating the cost and complexity associated with traditional product development and prototyping.

The Dimatix Materials Cartridge is a cartridge-based inkjet printhead used with the DMP. Based on FUJIFILM **Dimatix's proprietary Silicon MEMS** technology, the Dimatix Materials Cartridge is designed for highresolution, non-contact jetting of functional uids in a broad range of applications. The cartridge can deposit features to fabricate products such as organic thin- Im transistors (TFTs) and printed circuits. In biotechnology, the Dimatix Materials Cartridge allows researchers to closely pack large numbers of elements in DNA arrays, to permit more accurate and e cient analyses.

FUJIFILM Value from Innovation



Corporate Office: FUJIFILM Dimatix, Inc. 2250 Martin Avenue Santa Clara, CA 95050 USA

Tel: (408) 565-9150 Fax: (408) 565-9151 Email: info@dimatix.com New Hampshire Facility: FUJIFILM Dimatix, Inc. 109 Etna Road Lebanon, NH 03766 USA

Tel: (603) 443-5300 Fax: (603) 448-9870 m Email: info@dimatix.com Japan O ce: FUJIFILM Global Graphic Systems Co., Ltd. 2-26-30 Nishiazabu Minato-ku, Tokyo 106-0031 Japan Advanced Marking Strategy Division Phone: +81 3 6419 0530 Fax: +81 3 6419 9840 E-mail: dmp_ gs@ gs.fuji Im.co.jp Europe Office: Tel: +44 7739 863 505 Fax: +44 870 167 4328 Email: euro@dimatix.com

Korea Office:

Email: mdkorea@dimatix.com

Taiwan O ce: Email: mdtaiwan@dimatix.com China O ce: FUJIFILM Dimatix China Service Center Building 30, 1000 Jinhai Road Pudong New Area, Shanghai China 201206 Email: china@dimatix.com

Singapore O ce: Email: mdsingapore@dimatix.com

www.dimatix.com

PDS00085 Rev.05 3/11/19